



Email : sales@laboratoryinstrumentindia.com

Product Name :
Arsenator/Digital Arsenic Reader

Product Code :
LBNY-0028-220011



Description :

The displayed density measurement is correlated to a data table that corresponds directly to the concentration of arsenic in the water sample.

It allows the measure of arsenic concentrations ranging from 2 ppb to 100 ppb.

Portable, with a battery pack or used with an AC adaptor, the digital arsenic test scanner allows to simply, rapidly and accurately measure the colour density of the test strip used in the arsenic test kit.

Technical Specification :

Portable Arsenic test digital scanner to be used with visual arsenic test kit.



Laboratory Instrument India