



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Handpump,9-15m,Variant

**Product Code :**  
LBNY-0028-2200012



**Description :**

- 1 x Complete set of pump body with handle retainer and HDPE top bush
- 1 x Complete assembly of handle assembly with rubber cushions
- 1 x Complete assembly bottom flange assembly with bolts, with intermediate conical flange and suitable for bolting on to concrete platform
- 1 x Conical rubber grommet modified
- 1 x Double ended spanner M17 x M19
- 2 x Combination spanner
- 1 x MS retrieving rod, 2.20 meter long, electro galvanized
- 1 x Solvent cement 250mt in aluminium extruded bottle with PE stopper and sealed with plastic cap
- 4 x PVC pump rod with threaded joints cemented, 42mm OD x 3 meters
- 1 x PVC reducer (1 pc. PVC reducing socket 88.5mm ID and max. 93mm OD at the top side to fit with PVC casing pipe 89 mm OD x 60mm OD at the bottom side to fit with PVC cylinder reducer of 60mm ID & 1 pc. PVC cylinder reducer 60mm ID at the top side to fit with PVC reducing socket 60mm OD x 48mm ID at the bottom side to fit with 48 mm OD PVC lower casing)
- 1 x PVC top raising main with threaded joint cemented, 60mm OD x 3 meters
- 1 x PVC well screen assembly, 2 meters long
- 2 x PVC lower casing 48mm OD x 3 meters with bell connectors to fit item 15
- 1 x PVC sand trap complete assembly with end cap 48mm OD x 1 meter with ball connector to fit.
- 1 x PVC top pump rod with threaded joints cemented, 42mm OD x 3 meters
- 4 x PVC raising mains with threaded joints cemented, 60mm x 3 meters
- 1 x Complete assembly of 60mm OD x 2 meters long PVC cylinder pipe assembly with one complete foot valve assembly and one rubber seat.

---

**Technical Specification :**

Direct Action (low lift) hand pump, package Variant 2, for cylinder installation depths 9 to 12 meters.



**Laboratory Instrument India**