



Email : sales@laboratoryinstrumentindia.com

Product Name :
Electrode Dry Oven

Product Code :
LBNY-0004-820063



Description :

Electrode Dry Oven

Technical Specification :

Electrode Dry Oven

Temperature range: RT+10-250oC (Recommended temperature; 180oc)

Temperature error range: ± 1 oC

Temperature control unit: 0.1oC

Temperature uniformity: *35%

Inner chamber: stainless steel

Capacity 30L

Shelf 6Ea

Wattage: 0.8Kw, AC220V/3.6A/50-60Hz

Exhaust port: $\varnothing 28$ mm, top

Temperature control Mode: Dual Temperature segment Intelligent PID measured temperature, the set temperature and time are displayed at LCD screen, sensor: pt100



Laboratory Instrument India

