



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Water tank,collapsible,20000l,w/dist.kit

**Product Code :**

LBNY-0028-2200014



**Description :**

Coated with PVC suitable for drinking water storage.

Water tank, 20000 litres, collapsible, made from polyester

Complete set, includes 4 water distribution kits with 2 x 6 taps, connection kit, installation kit and repair kit. Food grade material.

**Technical Specification :**

Water tank, collapsible, 20 cbm (20,000 litres), made from polyester fibres and coated with PVC, UV proof, thermal stability -30°/+70°C, suitable for chlorinated drinking water storage (food grade material).

**Material:**

Mildew and rot resistant.

Ultra violet resistant.

Tear strength: 700 (W)/630 (F) N/5 cm

Tensile strength: 5075 (W)/5075 (F) N/5 cm

Adhesion: Bigger than 100 N/5 cm

Opaque - to inhibit algae growth.

Assembly: Heat welding/high frequency welding.

Weight: 930 g/m<sup>2</sup>.

Flame retardancy: None.

Temperature resistance: -30°C / +70°C.



**Laboratory Instrument India**