



Email : sales@laboratoryinstrumentindia.com

Product Name :
Handpump,SDWP,SS-Riser

Product Code :
LBNY-0028-2200015



Description :

Riser pipe, made of AISI 304 stainless steel, CrNi 18 10.

In length of 3 metres.

Nominal diameter: 32mm.

Nominal diameter 32mm (1"1/4), both ends male threaded, one end fitted with one female threaded coupling, the other end protected with a HDPE thread protector covering the entire length of the thread.

Each riser pipe comes complete with:

One Connecting Rod

Made of Stainless Steel, with M12 rolled male thread at one end, and M12 SS female threaded coupling welded at the other end.

One rod centralizer made of nitrile rubber or similar.

Nominal diameter: 32mm.

SS Rod Diameter: 10.8mm,

Length of 3 meters (+10/-0mm).

Technical Specification :

The Riser Pipe is supplied with:

1 x 3 meters SS connecting rod.

1 x centralizer riser pipe.

1 x hemp, 0.2 kg hank, for 100 riser pipes.

1 x tubes of alimentary grease for connecting rods.

Weight/Volume/Dimensions:

Estimated volume: 3.5cdm

Estimated weight: 11kg.



Laboratory Instrument India