



Email : sales@laboratoryinstrumentindia.com

Product Name :

Latrine kit, raised, f. hard or wet ground

Product Code :

LBNY-0028-2200016



Description :

Raised latrine kit for use in urban areas and on hard or flooded ground where impossible to dig a pit, composed of an above ground base with sludge tank, access ramp and elevated latrine slab, supporting its latrine superstructure.

Technical Specification :

Raised latrine kit, complete with and comprising of an above ground based support and substructure, latrine slab, superstructure, access ramp and an arrangement for collection of sludge.

The kit is easy in assembling to provide a stable and durable installation for a block of two (2) raised latrine stances/cubicles in first phase emergencies. Used for communal latrines, health centres, feeding centres, schools etc., it is composed of following main components:

- 1/- Superstructure with lockable door, roof/cover and natural ventilation/aeration
- 2/- Access ramp (steps with hand rail/support and platform).
- 3/- Substructure and excreta holding container/tank
- 4/- Squatting plate

This complete raised latrine cabin kit is for deployment in urban environment, hard/rocky ground or wet/flooded area installation where it is not possible or appropriate to dig a normal simple pit/trench.

The overall structures shall be of lightweight, rigid, durable and sturdy design, climatically stable and resistant to withstand heat, wind, dust, rain and manual cleaning, including stronger windy conditions.

The kit shall be complete, ready to install, and not requiring any materials to be sourced in country for the initial installation. (However it is allowed to foresee the design for later repairs with local available materials).



Laboratory Instrument India