



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Pump centrifugal,diesel,30cbm/hr@20m

**Product Code :**  
LBNY-0028-220021



**Description :**

Self-priming centrifugal pumping set, capacity 30 Cubic Metres per hour at around 2900 RPM, at 20 metre TMH (Total Manometric Head) for 1-m suction head.

Maximum suction lift 5 metres.

Centrifugal pump, diesel driven, with a capacity of 30 m<sup>3</sup>/hr at 20 m Total Manometric Head.

**Technical Specification :**

The pump set is supplied with following accessories:

1 x suction hose, semi-rigid, minimum 5 m.

1 x non return valve.

1 x tool kit for installation, operation and maintenance.

1 x delivery hose, minimum 20 m.

1 x stainless steel strainer.

1 x oil container, 4 litres, for diesel engine.

1 x spare parts kit for 2000 hours of normal operation, including various filters, seals and one injector.

Weight/Volume/Dimensions:

Estimated volume: 1.0 m<sup>3</sup>

Estimated weight: 120 kg.

---

# Laboratory Instrument India

