



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Air Compressor

**Product Code :**  
LBNY-0004-820064



**Description :**

Air Compressor

**Technical Specification :**

Air Compressor

Specifications

Capacity: 150 L/min(2.5HP)

Air production: (150-200)L/min

Weight/pressure: 21Kg/8bar

800W, 220V,50Hz

Accessory

Air Gun 2Ea

Air Hose(7.5m) 2Ea



**Laboratory Instrument India**