



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Bucket,HDPE,with lid and tap,14 l

**Product Code :**  
LBNY-0028-2200025



**Description :**

Bucket made of HDPE (High Density Poly Ethylene), tough durable UV resistant and safe for food and water storage.

The walls should meet the bottom of the bucket with a curved inside surface to prevent dirt accumulation and facilitate cleaning.

The bucket capacity is 14 litres.

Height 33cm, top diameter 31cm.

The top is reinforced to prevent ovaling.

**Technical Specification :**

The bucket has a tight-fitting lid of the same material of the bucket with an attached push-on cap.

Each bucket is supplied with a loose, sturdy plastic tap, stored during shipment in the bottom of each bucket.

The tap and its connection must be leak tight. At destination, the tap can be put in the predisposed hole in the bucket from the outside and fixed by screwing the locking nut from the inside.



**Laboratory Instrument India**

