



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Water tank,collapsible,10000l,w/dist.kit

**Product Code :**

LBNY-0028-2200028



**Description :**

Water tank, collapsible, 10 cbm (10,000 litres), made from polyester and coated with PVC, UV proof, thermal stability  $-30^{\circ}/+70^{\circ}\text{C}$ , suitable for chlorinated drinking water storage (food grade material).

Water tank, 10,000 litres, collapsible, suitable for drinking water storage, complete set including 2 water distribution kits with 2 x 6 taps, connection kit, installation kit and repair kit. Food grade material.

**Technical Specification :**

Material:

Mildew and rot resistant

Opaque - to inhibit algae growth

Ultra violet resistant

Assembly: Heat welding/high frequency welding.

Outer coating: Abrasion and UV-resistant fabric such as PVC, weight  $480\text{g}/\text{m}^2$ .

Inner coating: Black polyurethane, PVC or equivalent coating able to withstand, without alteration of the stored drinking water, chlorine concentration of  $0.5\text{ mg}/\text{litre}$  at a water temperature of  $35^{\circ}\text{C}$ , and peak chlorine concentration of  $5\text{ mg}/\text{litre}$  during cleaning process. Weight  $200\text{g}/\text{m}^2$ . Base fabric (substrate): Bi-composite polyester or equivalent material, weight  $250\text{ g}/\text{m}^2$ .

Threads per cm: 16 x 16.

Expected life time minimum three years as to chlorine and UV resistance.

Weight:  $930\text{ g}/\text{m}^2$ .

Adhesion: Bigger than 100 N/5 cm.

Flame retardant: None.

Tensile strength: 5075 (W)/5075 (F) N/5 cm.

---

Tear strength: 700 (W)/630 (F) N/5 cm.  
Temperature resistance: -30°C / +70°C.



**Laboratory Instrument India**