



Email : sales@laboratoryinstrumentindia.com

Product Name :
NaDCC,Chlorine Granules 56%, 1/2kg

Product Code :
LBNY-0028-2200030



www.laboratoryinstrumentindia.com

Description :

Chlorine, NaDCC 56%, for Water purification, as Sodium Dichloroisocyanurate (NaDCC) granules, in plastic pot of 0.5kg. Supplied with measuring spoon.

Chlorine as 100% Sodium Dichloroisocyanurate Dihydrate (NaDCC) in white granules, which releases 56% of available Chlorine:

Suitable for treatment of water for human consumption (drinking water purification).

Volume 20ml (= 15g of granules),

Supplied with measuring spoon, short handle,

Shelf life (expected): 3 years (Unopened container).

Technical Specification :

Chlorine, NaDCC 56%, for Water purification, as Sodium Dichloroisocyanurate (NaDCC) granules, in plastic pot of 0.5kg. Supplied with measuring spoon.



Laboratory Instrument India