



Email : sales@laboratoryinstrumentindia.com

Product Name :
Multipurpose Cloth, cotton, 1x1.5m

Product Code :
LBNY-0028-220034



Description :

Multipurpose Cloth.

Wrap: minimum 21 threads per cm,

Weft: minimum 18.5 threads per cm.

Color: natural, light brown or grey.

Cloth made of 100% cotton with approximate size of 1 m x 1.5 m and fabric weight of 170-180g/m².

Technical Specification :

Multipurpose Cloth, cotton, 1m x 1.5m.



Laboratory Instrument India