



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Multipurpose Cloth, cotton, 1x1.5m

**Product Code :**  
LBNY-0028-220034



**Description :**

Multipurpose Cloth.

Wrap: minimum 21 threads per cm,

Weft: minimum 18.5 threads per cm.

Color: natural, light brown or grey.

Cloth made of 100% cotton with approximate size of 1 m x 1.5 m and fabric weight of 170-180g/m<sup>2</sup>.

**Technical Specification :**

Multipurpose Cloth, cotton, 1m x 1.5m.



**Laboratory Instrument India**