



Email : sales@laboratoryinstrumentindia.com

Product Name :
Water-cont,LDPE,10l,collapsibl.

Product Code :
LBNY-0028-2200035



Description :

Collapsible water container with a capacity of 10 litres and a 50mm opening, made from food grade material. Without logo.

Technical Specification :

Material

It is designed for minimum 6 months of use under tropical conditions. The collapsible water container is made of virgin food-grade plastic (e.g., Low-Density Polyethylene, High-Density Polyethylene or Polypropylene), UV resistant and safe for food and water storage.

Design

The 10L collapsible water container comes with a robust and comfortable handle on one side. The opening has a diameter of at least 50mm, with a matching screw cap attached with a retaining strap. The collapsible 10 L water container is easy to fill up and pour, comfortable to carry for the end user and comes with no sharp edges. It is semi-transparent to identify the level of liquid inside visually.

Impact Resistance / Drop Test

The collapsible jerrycan resists a minimum of 3 consecutive drops from 2m height, dropped from different angles.

Laboratory Instrument India

