



Email : sales@laboratoryinstrumentindia.com

Product Name :
Bucket,HDPE,with lid,14 l

Product Code :
LBNY-0028-2200036



Description :

Bucket made of HDPE (High Density Poly Ethylene), UV resistant and safe for food and water storage. The walls shall meet the bottom of the bucket with a curved surface to prevent dirt accumulation and facilitate cleaning.

The bucket has a tight-fitting lid of the same material of the bucket with an attached push-on cap.

The bucket capacity is 14 litres.

Height 33 cm, top diameter 31 cm

The top is reinforced to prevent ovaling.

Technical Specification :

Weight/Volume/Dimensions:

Estimated volume: 15cdm

Estimated dimensions: Top diameter 31cm, height 33cm

Estimated weight: 1.1kg



Laboratory Instrument India

