



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Bucket,HDPE, w/lid,10cm cap and tap,14 l

**Product Code :**

LBNY-0028-2200037



**Description :**

Bucket designed solely for transport and storage of drinking water, made of food grade virgin HDPE (High Density PolyEthylene)/virgin LDPE (Low Density PolyEthylene), tough durable UV-resistant that will not deteriorate in harsh sunlight, and safe for food and drinking water storage.

Bucket (container), stackable, round, made of HDPE (High Density Poly Ethylene), with lid, large cap, large plastic handle and tap, capacity 14 litres.

**Technical Specification :**

The bucket capacity is 14 litres min.

Top diameter 300/310 mm  $\pm$ 5%,

Bottom diameter 240 mm  $\pm$ 5%.

Height 300 mm,

Indicated dimensions are approximate.



**Laboratory Instrument India**