



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Bucket, HDPE, w/lid and 10cm cap, 14 l

**Product Code :**

LBNY-0028-2200038



**Description :**

Bucket designed solely for transport and storage of drinking water, made of food grade virgin HDPE (High Density PolyEthylene)/virgin LDPE (Low Density PolyEthylene), tough durable UV-resistant that will not deteriorate in harsh sunlight, and safe for food and drinking water storage.

**Technical Specification :**

The bucket capacity is 14 litres min.  
Top diameter 300/310 mm  $\pm 5\%$ ,  
Bottom diameter 240 mm  $\pm 5\%$ .  
Height 300 mm.



**Laboratory Instrument India**