



Email : sales@laboratoryinstrumentindia.com

Product Name :

Bucket, HDPE, w/lid and 10cm cap, 14 l

Product Code :

LBNY-0028-2200038



Description :

Bucket designed solely for transport and storage of drinking water, made of food grade virgin HDPE (High Density PolyEthylene)/virgin LDPE (Low Density PolyEthylene), tough durable UV-resistant that will not deteriorate in harsh sunlight, and safe for food and drinking water storage.

Technical Specification :

The bucket capacity is 14 litres min.
Top diameter 300/310 mm $\pm 5\%$,
Bottom diameter 240 mm $\pm 5\%$.
Height 300 mm.



Laboratory Instrument India