



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pan and "P" trap, for sq. plate

Product Code :
LBNY-0028-220041



Description :

The pan is moulded in PP and can be bolted to the underside of squatting plate using the fasteners provided. The trap is fixed to the pan using the adhesive that has been provided. There are 2 types of traps, a "P" trap and a Goose neck trap, both of which are moulded in HDPE.

Technical Specification :

Weight/Volume/Dimensions:

For pan:

Estimated volume: 40 cdm

Estimated dimensions: 470 x 250 x 345 mm

Estimated weight: 0.5 kg.



Laboratory Instrument India