



Email : sales@laboratoryinstrumentindia.com

Product Name :
Footpump,80m well

Product Code :
LBNY-0028-2200043



Description :

Complete Vergnet hydro-pump; hydraulic transmission, foot operated and multi-operator pump type, supplied with 80 metres each of complete pressure (command/transmission) and riser (discharge) pipes.

Technical Specification :

Supplied with
1 x pump barrel
1 x lower pumping element, stainless steel
1 x delivery hose
1 x pack of clips
1 x installation manual
1 x base plate with sealing pins
1 x command cylinder, stainless steel
1 x command hose
1 x maintenance manual
1 x tool kit



Laboratory Instrument India

