



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Footpump,60m well

**Product Code :**  
LBNY-0028-2200046



**Description :**

Complete Vergnet hydro-pump; hydraulic transmission, foot operated and multi-operator pump type, supplied with 60 metres each of complete pressure (command/transmission) and riser (discharge) pipes.

Supplied with

- 1 x pump barrel
- 1 x lower pumping element, stainless steel
- 1 x delivery hose
- 1 x pack of clips
- 1 x base plate with sealing pins
- 1 x command cylinder, stainless steel
- 1 x command hose
- 1 x maintenance manual
- 1 x installation manual
- 1 x tool kit.

**Technical Specification :**

Weight and Volume  
Estimated Volume: 0.8m<sup>3</sup>  
Estimated Weight: 90kg

---

# Laboratory Instrument India

