



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Double Open End Spanner Chrome Vanadium

**Product Code :**  
LBNY-0029-220001



**Description :**

Drop forged from high quality Chrome Vanadium Steel  
Available in Chrome Mirror Finish, Satin Finish & Matt Finish.  
Induction Heat treated.

**Technical Specification :**

<b>SIZE (MM)</b>	<b>NET WT. IN GMS.</b>	<b>LENGTH (MM)</b>	<b>QTY. / CTN. (PCS.)</b>	<b>NET WT. / CTN.(KGS)</b>	<b>GR.</b>
6x7	24	122	1000	24	25.5
8x9	32	140	600	19.2	20.7
10x11	55	155	400	22	23.5
12x13	76	171	300	22.8	24.3
14x15	112	190	200	22.4	23.9
16x17	130	204	150	19.5	21
18x19	165	220	100	16.5	18
20x22	210	235	100	21	22.5
21x23	240	248	80	19.2	20.7
24x26	315	267	60	18.9	20.4
24x27	315	267	60	18.9	20.4
25x28	360	280	60	21.6	23.1
30x32	485	301	40	19.4	20.9



**Laboratory Instrument India**