



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Bi-Hexagonal Ring Spanner Elliptical Type

**Product Code :**  
LBNY-0029-220002



**Description :**

Drop forged from high quality Chrome Vanadium Steel  
Available in Chrome Mirror Finish, Satin Finish, Matt Finish.  
Induction Heat treated

**Technical Specification :**

SIZE (MM)	NET WT. IN GMS.		LENGTH (MM)		QTY. / CTN. (PCS.)	NET WT. / CTN.(KGS)	GR.
6x7	50	183	400	20	21.5		
8x9	70	192	300	21	22.5		
10x11	92	202	250	23	24.5		
12x13	125	217	200	25	26.5		
14x15	150	233	150	22.5	24		
16x17	170	251	100	17	18.5		
18x19	210	268	100	21	22.5		
20x22	305	296	60	18.3	19.8		
21x23	325	322	60	19.5	21		
24x27	395	330	40	15.8	17.3		
25x28	425	345	40	17	18.5		
30x32	560	370	30	16.8	18.3		



**Laboratory Instrument India**