



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
ELLIPTICAL COMBINATION SPANNERS

**Product Code :**  
LBNY-0029-220006



**Description :**

Drop forged from high quality Chrome Vanadium Steel  
Available in Chrome Mirror Finish, Satin Finish & Matt Finish.  
Induction Heat treated

**Technical Specification :**

SIZE (MM)	SIZE (INCH)	NET WT. IN GMS.	LENGTH (MM)	QTY. / CTN. (PCS.)	NET WT. / CTN.(KGS)
6	1/4	18	100	1000	18
7	--	22	110	1000	22
8	5/6	26	120	1000	26
9	--	30	130	800	24
10	3/8	42	140	600	25.2
11	7/16	52	150	400	20.8
12	--	60	160	400	24
13	1/2	72	170	300	21.6
14	9/16	86	180	250	21.5
15	--	92	190	200	18.4
16	5/8	109	200	200	21.8
17	11/16	132	210	200	26.4
18	--	152	220	150	22.8
19	3/4	160	230	125	20
20	--	182	240	125	22.75

---

21	13/16	192	250	100	19.2	20.7
22	7/8	210	260	100	21	22.5
23	--	250	270	80	20	21.5
24	15/16	260	280	80	20.8	22.3
25	1"	300	290	60	18	19.5
26	--	320	300	60	19.2	20.7
27	1.1/16	390	310	60	23.4	24.9
28	1.1/8	435	320	40	17.4	18.9
29	--	435	330	40	17.4	18.9
30	--	460	340	40	18.4	19.9
32	1.1/4	560	360	40	22.4	23.9



**Laboratory Instrument India**