



Email : sales@laboratoryinstrumentindia.com

Product Name :
Staff Gauges Plates + Accessories

Product Code :
LBNY-0002-830003



Description :

Staff Gauges Plates + Accessories

Technical Specification :

Staff Gauge plates are primarily used to undertake manual readings of water level in water bodies. The plates are made from corrosion resistant aluminum with a stove baked enamel background. They are supplied in screen printed markings, with graduations of one centimetre. The plates can be provided with or without a range of metre graduations (see images). They are available flat or in a rounded shape to allow fixing to 80 mm metal pipe posts.



Laboratory Instrument India