



Email : sales@laboratoryinstrumentindia.com

Product Name :
DOUBLE OPEN ENDED SPANNERS (COLD STAMPED)

Product Code :
LBNY-0029-2200010



Description :

Drop forged from high quality Chrome Vanadium Steel
Double Layer Chrome Plated / Satin Finish.
Anti Slip panel, firm grip even with oily hands.

Technical Specification :

SIZE (MM)	NET WT. IN GMS.	LENGTH (MM)	QTY. / CTN. (PCS.)	NET WT. / CTN.(KGS)	GR.
6x7	23	122	1000	23.00	24.50
8x9	32	140	600	19.20	20.70
10x11	50	155	400	20.00	21.50
12x13	62	171	300	18.60	20.10
14x15	80	190	200	16.00	17.50
16x17	110	204	150	16.50	18.00
18x19	150	220	100	15.00	16.50
20x22	200	235	100	20.00	21.50
21x23	230	248	80	18.40	19.90
24x26	320	267	60	19.20	20.70
24x27	275	267	60	16.50	18.00
25x28	345	280	50	17.25	18.75
30x32	450	301	40	18.00	19.50



Laboratory Instrument India