



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
RING SPANNERS (COLD STAMPED)

**Product Code :**  
LBNY-0029-2200011



**Description :**

Drop forged from high quality Chrome Vanadium Steel  
Double Layer Chrome Plated / Satin Finish.  
Anti Slip panel, firm grip even with oily hands.

**Technical Specification :**

SIZE (MM)	NET WT. IN GMS.	LENGTH (MM)	QTY. / CTN. (PCS.)	NET WT. / CTN.(KGS)	GR.
6x7	50	183	400	20.00	21.50
8x9	70	192	300	21.00	22.50
10x11	92	202	250	23.00	24.50
12x13	125	217	200	25.00	26.50
14x15	150	233	150	22.50	24.00
16x17	170	251	100	17.00	18.50
18x19	210	268	80	16.80	18.30
20x22	305	296	60	18.30	19.80
21x23	325	322	60	19.50	21.00
24x27	395	330	40	15.80	17.30
25x28	425	345	40	17.00	18.50
30x32	560	370	30	16.80	18.30



**Laboratory Instrument India**