



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
COMBINATION SPANNERS (COLD STAMPED)

**Product Code :**  
LBNY-0029-2200012



**Description :**

Drop forged from high quality Chrome Vanadium Steel  
Anti Slip panel, firm grip even with oily hands.

**Technical Specification :**

SIZE (MM)	SIZE(INCH)	NET WT. IN GMS.	LENGTH (MM)	QTY. / CTN. (PCS.)	NET WT. / CTN.(KGS)
6	1/4	18	100	1000	18.00 19.50
7	--	22	110	1000	22.00 23.50
8	5/6	26	120	800	20.80 22.30
9	--	30	130	600	18.00 19.50
10	3/8	42	140	500	21.00 22.50
11	7/16	52	150	400	20.80 22.30
12	--	60	160	300	18.00 19.50
13	1/2	72	170	200	14.40 15.90
14	9/16	86	180	200	17.20 18.70
15	--	92	190	150	13.80 15.30
16	5/8	109	200	100	10.90 12.40
17	11/16	132	210	100	13.20 14.70
18	--	152	220	100	15.20 16.70
19	3/4	160	230	80	12.80 14.30
20	--	182	240	60	10.92 12.42
21	13/16	192	250	60	11.52 13.02

---

22	7/8	210	260	60	12.60	14.10
23	--	250	270	60	15.00	16.50
24	15/16	260	280	50	13.00	14.50
25	1	300	290	40	12.00	13.50
26	--	320	300	40	12.80	14.30
27	1.1/16	390	310	30	11.70	13.20
28	1.1/8	435	320	30	13.05	14.55
29	--	435	330	30	13.05	14.55
30	--	460	340	20	9.20	10.70
32	1.1/4	560	360	20	11.20	12.70



**Laboratory Instrument India**