



Email : sales@laboratoryinstrumentindia.com

Product Name :
DOUBLE OPEN END SPANNERS (RECESSED PANEL
DROP FORGED)

Product Code :
LBNY-0029-2200013



Description :

Anti Slip panel, firm grip even with oily hands
Drop forged from high quality Chrome Vanadium Steel, Carbon Steel

Technical Specification :

SIZE (MM)	NET WT. IN GMS.	LENGTH (MM)	NET WT. / CTN.(KGS)	GR. WT.(KGS)
6x7	20 98 20	21.5		
8x9	25 102 25	26.5		
10x11	35 110 28	29.5		
12x13	50 128 25	26.5		
14x15	65 150 19.5	21		
16x17	75 162 18.75	20.25		
18x19	110 181 22	23.5		
20x21	155 199 23.25	24.75		
22x23	165 208 16.5	18		
24x26	245 225 19.6	21.1		
24x27	245 225 19.6	21.1		
25x28	265 233 21.2	22.7		
30x32	325 248 19.5	21		



Laboratory Instrument India