



Email : sales@laboratoryinstrumentindia.com

Product Name :
DOUBLE OPEN END SPANNERS (RECESSED PANEL
DROP FORGED)

Product Code :
LBNY-0029-2200013



Description :

Anti Slip panel, firm grip even with oily hands
Drop forged from high quality Chrome Vanadium Steel, Carbon Steel

Technical Specification :

SIZE (MM)	NET WT. IN GMS.		LENGTH (MM)		NET WT. / CTN.(KGS)		GR. WT.(KGS)	
6x7	20	98	20	21.5				
8x9	25	102	25	26.5				
10x11	35	110	28	29.5				
12x13	50	128	25	26.5				
14x15	65	150	19.5	21				
16x17	75	162	18.75	20.25				
18x19	110	181	22	23.5				
20x21	155	199	23.25	24.75				
22x23	165	208	16.5	18				
24x26	245	225	19.6	21.1				
24x27	245	225	19.6	21.1				
25x28	265	233	21.2	22.7				
30x32	325	248	19.5	21				



Laboratory Instrument India