



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
COMBINATION SPANNERS (RECESSED PANEL DROP FORGED)

**Product Code :**  
LBNY-0029-2200015



**Description :**

Anti Slip panel, firm grip even with oily hands  
Drop forged from high quality Chrome Vanadium Steel, Carbon Steel.

**Technical Specification :**

SIZE (MM)	SIZE(INCH)	NET WT. IN GMS.	LENGTH (MM)	NET WT. / CTN.(KGS)	GR. WT.(KGS)
6	1/4	18	100	21.6	23.1
7	--	20	110	20	21.5
8	5/6	24	120	19.2	20.7
9	--	28	130	800	22.4 23.9
10	3/8	38	140	600	22.8 24.3
11	7/16	48	150	400	19.2 20.7
12	--	55	160	400	22 23.5
13	1/2	60	170	350	21 22.5
14	9/16	80	180	300	24 25.5
15	--	90	190	250	22.5 24
16	5/8	105	200	200	21 22.5
17	11/16	115	210	200	23 24.5
18	--	130	220	150	19.5 21
19	3/4	145	230	125	18.125 19.625
20	--	170	240	125	21.25 22.75
21	13/16	175	250	125	21.875 23.375

---

22	7/8	185	260	125	23.125	24.625
23	--	230	260	100	23	24.5
24	15/16	260	270	80	20.8	22.3
25	1"	300	280	60	18	19.5
26	--	325	280	60	19.5	21
27	1.1/16	350	290	60	21	22.5
28	1.1/8	400	300	60	24	25.5
29	--	455	300	40	18.2	19.7
30	--	480	310	40	19.2	20.7
32	1.1/4	575	320	40	23	24.5



**Laboratory Instrument India**