



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
RING SPANNERS (RECESSED PANEL CHROME  
VANADIUM)

**Product Code :**  
LBNY-0029-2200017



**Description :**

Drop forged from high quality Chrome Vanadium Steel, Carbon Steel  
Anti Slip panel, firm grip even with oily hands.

**Technical Specification :**

<b>SIZE (MM)</b>	<b>NET WT. IN GMS.</b>		<b>LENGTH (MM)</b>		<b>NET WT. / CTN.(KGS)</b>		<b>GR. WT.(KGS)</b>	
6x7	50	183	20	21.5				
8x9	60	192	18	19.5				
10x11	85	202	21.25	22.75				
12x13	115	217	23	24.5				
14x15	130	233	19.5	21				
16x17	170	251	17	18.5				
18x19	200	268	20	21.5				
20x22	275	296	20.625	22.125				
21x23	295	322	17.7	19.2				
24x27	370	330	18.5	20				
25x28	395	345	15.8	17.3				
30x32	515	370	10.3	11.8				

---

# Laboratory Instrument India

