



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
COMBINATION SPANNERS (RAISED PANEL)

**Product Code :**  
LBNY-0029-2200019



**Description :**

Double Layer Chrome Plated / Satin Finish

**Technical Specification :**

SIZE (MM)	SIZE(INCH)	NET WT. IN GMS.	LENGTH (MM)	NET WT. / CTN.(KGS)	GR. WT.(KGS)
6	1/4	16	95	22.4	23.9
7	--	20	97	24	25.5
8	5/6	28	108	22.4	23.9
9	--	30	114	24	25.5
10	3/8	40	128	24	25.5
11	7/16	43	132	21.5	23
12	--	54	137	21.6	23.1
13	1/2	65	145	19.5	21
14	9/16	80	164	24	25.5
15	--	90	175	22.5	24
16	5/8	100	185	20	21.5
17	11/16	135	205	20.25	21.75
18	--	150	205	22.5	24
19	3/4	160	232	20	21.5
20	--	170	232	21.25	22.75
21	13/16	190	258	19	20.5
22	7/8	215	289	21.5	23
23	--	225	289	22.5	24

---

24	15/16	260	313	20.8	22.3
25	1"	270	324	21.6	23.1
26	--	300	324	24	25.5
27	1.1/16	350	361	21	22.5
28	1.1/8	400	369	24	25.5
30	--	525	381	21	22.5
32	1.1/4	620	391	24.8	26.3



**Laboratory Instrument India**