



Email : sales@laboratoryinstrumentindia.com

Product Name :
C-TYPE HALF MOON SPANNERS RING

Product Code :
LBNY-0029-2200020



Description :

These are drop forged and made of high quality chrome vanadium steel.
These spanners are widely used for its durability and quality.
Chrome Vanadium Steel, Full Polished or Satin Finish

Technical Specification :

SIZE	LENGTH	NET WEIGHT	NET WT. / CTN.(KGS)	GR. WT.(KGS)
8X10	115	40	20	21.5
10X11	127	50	20	21.5
11x13	150	75	22.5	24
12X13	150	73	18.25	19.75
13x15	165	105	21	22.5
14X15	165	108	21.6	23.1
13x17	185	110	22	23.5
14x17	185	120	24	25.5
15X17	185	120	18	19.5
16X17	185	125	18.75	20.25
17X19	185	125	18.75	20.25
18X19	202	150	18.75	20.25
19x21	220	200	20	21.5
19x22	220	200	20	21.5
20X22	220	205	20.5	22
21X23	230	205	20.5	22



Laboratory Instrument India