



Email : sales@laboratoryinstrumentindia.com

Product Name :
C-TYPE HALF MOON SPANNERS RING

Product Code :
LBNY-0029-2200020



Description :

These are drop forged and made of high quality chrome vanadium steel.
These spanners are widely used for its durability and quality.
Chrome Vanadium Steel, Full Polished or Satin Finish

Technical Specification :

SIZE	LENGTH	NET WEIGHT	NET WT. / CTN.(KGS)	GR. WT.(KGS)
8X10	115 40	20 21.5		
10X11	127 50	20 21.5		
11x13	150 75	22.5 24		
12X13	150 73	18.25 19.75		
13x15	165 105	21 22.5		
14X15	165 108	21.6 23.1		
13x17	185 110	22 23.5		
14x17	185 120	24 25.5		
15X17	185 120	18 19.5		
16X17	185 125	18.75 20.25		
17X19	185 125	18.75 20.25		
18X19	202 150	18.75 20.25		
19x21	220 200	20 21.5		
19x22	220 200	20 21.5		
20X22	220 205	20.5 22		
21X23	230 205	20.5 22		



Laboratory Instrument India