



Email : sales@laboratoryinstrumentindia.com

Product Name :
C-TYPE HALF MOON SPANNERS COMBINATION

Product Code :
LBNY-0029-2200021



Description :

These are drop forged and made of high quality chrome vanadium steel.
These spanners are widely used for its durability and quality.
Chrome Vanadium Steel, Full Polished or Satin Finish.

Technical Specification :

SIZE	SAE	LENGTH	NET WEIGHT	NET WT. / CTN.(KGS)	GR. WT.(KGS)
10	3/8"	140	55	22	23.7
11	7/16"	140	51	20.4	22.1
12	**	155	84	16.8	17.9
13	1/2"	155	83	16.6	17.7
14	9/16"	170	100	20	21.1
15	**	170	95	19	20
16	5/8"	185	130	19.5	20.5
17	11/16"	185	130	19.5	20.9
19	3/4"	205	170	17	18.1
22		19.5	20.4		

Laboratory Instrument India

