



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
S-TYPE RING SPANNERS

**Product Code :**  
LBNY-0029-220022



**Description :**

Optimum allow-steel structure and head treatment throughout ensure high strength. the long shank of wrench allows extra leverage without compromising safety.

**Technical Specification :**

SIZE(MM)	LENGTH(MM)	NET WT. / CTN.(KGS)	GR. WT.(KGS)
8X9	20 22		
10X11	128 16.8	18.2	
12X13	142 21.6	23	
14X15	165 17.25	18.2	
16X17	195 16.5	17.6	
18X19	222 17.1	18.2	
21X22	20.475	21.425	



**Laboratory Instrument India**