



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
S-TYPE COMBINATION SPANNERS

**Product Code :**  
LBNY-0029-2200023



**Description :**

Optimum allow-steel structure and head treatment throughout ensure high strength. the long shank of wrench allows extra leverage without compromising safety.

**Technical Specification :**

SIZE(MM)	SAE	LENGTH(MM)	NET WT.IN GMS.	NET WT. / CTN.(KGS)	GR. WT. (KGS)
8	20.5	22.5	0.025		
10	3/8"	135	19.8	21.2	0.024
11	7/16"	135	18	19.4	0.024
12	**	150	17.4	18.5	0.022
13	1/2"	150	17.8	18.9	0.022
14	9/16"	175	20.25	21.2	0.019
15	**	175	19.2	20.05	0.019
16	5/8"	205	19.5	20.5	0.022
17	11/16"	205	19.5	20.6	0.022
19	3/4"	230	19.125	20.075	0.02



**Laboratory Instrument India**

