



Email : sales@laboratoryinstrumentindia.com

Product Name :
TORX SCREW WRENCH

Product Code :
LBNY-0029-2200024



Description :

Double Ended Ring Wrenches Straight for outside Torx Screws.

Technical Specification :

SIZE(MM)	NET WT.IN GMS.	LENGTH(MM)	NET WT. / CTN.(KGS)	GR. WT.(KGS)
6X7	12 95	21.6	23.1	
8X9	20 110	20	21.5	
10X11	30 125	21	22.5	
12X13	50 140	20	21.5	
14X15	60 156	24	25.5	
16X17	60 170	24	25.5	
17X19	95 180	23.75	25.25	
18X19	90 180	22.5	24	
20X22	138 200	20.7	22.2	
21X23	157 220	23.55	25.05	
24X26	290 272	23.2	24.7	
24X27	290 272	23.2	24.7	
25X28	285 270	22.8	24.3	
30X32	405 305	20.25	21.75	

Laboratory Instrument India

