



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
SINGLE OPEN END SPANNERS

**Product Code :**  
LBNY-0029-2200026



**Description :**

Fully hardened & tempered  
Drop forged from carbon steel  
Chrome plated/zinc plated/black phosphate finish.

**Technical Specification :**

| SIZE(MM) | NET WT.IN GMS. | LENGTH(MM) | NET WT. / CTN.(KGS) | GR. WT.(KGS) |
|----------|----------------|------------|---------------------|--------------|
| 6        | 10             | 75         | 20.00               | 21.50        |
| 7        | 11             | 78         | 22.00               | 23.50        |
| 8        | 16             | 95         | 22.40               | 23.90        |
| 9        | 23             | 105        | 20.70               | 22.20        |
| 10       | 23             | 105        | 20.70               | 22.20        |
| 11       | 23             | 111        | 20.70               | 22.20        |
| 12       | 32             | 125        | 22.40               | 23.90        |
| 13       | 35             | 125        | 21.00               | 22.50        |
| 14       | 46             | 140        | 23.00               | 24.50        |
| 15       | 50             | 142        | 25.00               | 26.50        |
| 16       | 55             | 157        | 22.00               | 23.50        |
| 17       | 65             | 160        | 19.50               | 21.00        |
| 18       | 65             | 160        | 19.50               | 21.00        |
| 19       | 82             | 170        | 24.60               | 26.10        |
| 20       | 80             | 170        | 24.00               | 25.50        |
| 21       | 120            | 195        | 18.00               | 19.50        |

---

|     |      |      |       |       |
|-----|------|------|-------|-------|
| 22  | 125  | 195  | 18.75 | 20.25 |
| 24  | 155  | 215  | 23.25 | 24.75 |
| 25  | 155  | 215  | 23.25 | 24.75 |
| 27  | 200  | 240  | 24.00 | 25.50 |
| 30  | 270  | 260  | 21.60 | 23.10 |
| 32  | 305  | 275  | 24.40 | 25.90 |
| 34  | 310  | 275  | 24.80 | 26.30 |
| 36  | 425  | 300  | 21.25 | 22.75 |
| 38  | 445  | 305  | 22.25 | 23.75 |
| 41  | 615  | 340  | 24.60 | 26.10 |
| 46  | 755  | 375  | 22.65 | 24.15 |
| 50  | 1085 | 405  | 21.70 | 23.20 |
| 55  | 1280 | 450  | 20.48 | 21.98 |
| 60  | 1522 | 475  | 18.26 | 19.76 |
| 65  | 1800 | 515  | 21.60 | 23.10 |
| 70  | 2205 | 560  | 22.05 | 23.55 |
| 75  | 3180 | 600  | 25.44 | 26.94 |
| 80  | 3930 | 670  | 23.58 | 25.08 |
| 85  | 3980 | 680  | 23.88 | 25.38 |
| 90  | 3700 | 680  | 22.20 | 23.70 |
| 95  | 5800 | 847  | 23.20 | 24.70 |
| 100 | 5500 | 847  | 22.00 | 23.50 |
| 105 | 7450 | 1000 | 14.90 | 16.40 |
| 110 | 7400 | 1000 | 14.80 | 16.30 |
| 115 | 7300 | 1000 | 14.60 | 16.10 |
| 120 | 7250 | 1000 | 14.50 | 16.00 |
| 125 | 7200 | 1000 | 14.40 | 15.90 |
| 130 | 7150 | 1000 | 14.30 | 15.80 |
| 135 | 7100 | 1000 | 14.20 | 15.70 |



**Laboratory Instrument India**