



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
SINGLE OPEN END SPANNERS

**Product Code :**  
LBNY-0029-2200026



**Description :**

Fully hardened & tempered  
Drop forged from carbon steel  
Chrome plated/zinc plated/black phosphate finish.

**Technical Specification :**

SIZE(MM)	NET WT.IN GMS.	LENGTH(MM)	NET WT. / CTN.(KGS)	GR. WT.(KGS)
6	10	75	20.00	21.50
7	11	78	22.00	23.50
8	16	95	22.40	23.90
9	23	105	20.70	22.20
10	23	105	20.70	22.20
11	23	111	20.70	22.20
12	32	125	22.40	23.90
13	35	125	21.00	22.50
14	46	140	23.00	24.50
15	50	142	25.00	26.50
16	55	157	22.00	23.50
17	65	160	19.50	21.00
18	65	160	19.50	21.00
19	82	170	24.60	26.10
20	80	170	24.00	25.50
21	120	195	18.00	19.50

---

22	125	195	18.75	20.25
24	155	215	23.25	24.75
25	155	215	23.25	24.75
27	200	240	24.00	25.50
30	270	260	21.60	23.10
32	305	275	24.40	25.90
34	310	275	24.80	26.30
36	425	300	21.25	22.75
38	445	305	22.25	23.75
41	615	340	24.60	26.10
46	755	375	22.65	24.15
50	1085	405	21.70	23.20
55	1280	450	20.48	21.98
60	1522	475	18.26	19.76
65	1800	515	21.60	23.10
70	2205	560	22.05	23.55
75	3180	600	25.44	26.94
80	3930	670	23.58	25.08
85	3980	680	23.88	25.38
90	3700	680	22.20	23.70
95	5800	847	23.20	24.70
100	5500	847	22.00	23.50
105	7450	1000	14.90	16.40
110	7400	1000	14.80	16.30
115	7300	1000	14.60	16.10
120	7250	1000	14.50	16.00
125	7200	1000	14.40	15.90
130	7150	1000	14.30	15.80
135	7100	1000	14.20	15.70



**Laboratory Instrument India**