



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
SLOGGING WRENCH OPEN SPANNERS

**Product Code :**  
LBNY-0029-2200029



**Description :**

Forged from high grade Selected Alloy with head set at 15 degree  
Heavy duty impact spanner with strong profiled handle and durable striking end.

**Technical Specification :**

<b>SIZE(MM)</b>	<b>NET WT.IN GMS.</b>	<b>LENGTH(MM)</b>	<b>NET WT. / CTN.(KGS)</b>	<b>GR. WT.(KGS)</b>
24	310	160	24.80	26.30
27	41	180	2.46	3.96
30	530	190	21.20	22.70
32	615	196	22.14	23.64
33	61	196	2.20	3.70
34	687	200	20.61	22.11
35	687	200	20.61	22.11
36	720	210	21.60	23.10
41	900	225	21.60	23.10
42	1090	230	21.80	23.30
46	1420	251	22.72	24.22
50	1600	275	19.20	20.70
55	2429	304	24.29	25.79
60	2800	320	22.40	23.90
65	3250	340	19.50	21.00
70	4000	365	24.00	25.50
75	5073	378	20.29	21.79

---

80	7000	400	28.00	29.50
85	6800	400	27.20	28.70
90	9600	445	19.20	20.70
95	9400	445	18.80	20.30
100	11020	485	22.04	23.54
105	11000	485	22.00	23.50
110	12700	500	25.40	26.90



## Laboratory Instrument India