



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
SLOGGING WRENCH OPEN SPANNERS

**Product Code :**  
LBNY-0029-2200029



**Description :**

Forged from high grade Selected Alloy with head set at 15 degree  
Heavy duty impact spanner with strong profiled handle and durable striking end.

**Technical Specification :**

| <b>SIZE(MM)</b> | <b>NET WT.IN GMS.</b> | <b>LENGTH(MM)</b> | <b>NET WT. / CTN.(KGS)</b> | <b>GR. WT.(KGS)</b> |
|-----------------|-----------------------|-------------------|----------------------------|---------------------|
| 24              | 310                   | 160               | 24.80                      | 26.30               |
| 27              | 41                    | 180               | 2.46                       | 3.96                |
| 30              | 530                   | 190               | 21.20                      | 22.70               |
| 32              | 615                   | 196               | 22.14                      | 23.64               |
| 33              | 61                    | 196               | 2.20                       | 3.70                |
| 34              | 687                   | 200               | 20.61                      | 22.11               |
| 35              | 687                   | 200               | 20.61                      | 22.11               |
| 36              | 720                   | 210               | 21.60                      | 23.10               |
| 41              | 900                   | 225               | 21.60                      | 23.10               |
| 42              | 1090                  | 230               | 21.80                      | 23.30               |
| 46              | 1420                  | 251               | 22.72                      | 24.22               |
| 50              | 1600                  | 275               | 19.20                      | 20.70               |
| 55              | 2429                  | 304               | 24.29                      | 25.79               |
| 60              | 2800                  | 320               | 22.40                      | 23.90               |
| 65              | 3250                  | 340               | 19.50                      | 21.00               |
| 70              | 4000                  | 365               | 24.00                      | 25.50               |
| 75              | 5073                  | 378               | 20.29                      | 21.79               |

---

|     |       |     |       |       |
|-----|-------|-----|-------|-------|
| 80  | 7000  | 400 | 28.00 | 29.50 |
| 85  | 6800  | 400 | 27.20 | 28.70 |
| 90  | 9600  | 445 | 19.20 | 20.70 |
| 95  | 9400  | 445 | 18.80 | 20.30 |
| 100 | 11020 | 485 | 22.04 | 23.54 |
| 105 | 11000 | 485 | 22.00 | 23.50 |
| 110 | 12700 | 500 | 25.40 | 26.90 |



## Laboratory Instrument India