



Email : sales@laboratoryinstrumentindia.com

Product Name :
GAS SPANNER - Z TYPE

Product Code :
LBNY-0029-2200031



Description :

Hardened and Tempered
Drop Forged from CARBON STEEL
Black or Chrome Plated.

Technical Specification :

SIZE(MM)	NET WT.IN GMS.	LENGTH(MM)	NET WT. / CTN.(KGS)	GR. WT.(KGS)
25 105 120	21	22.5		
27 125 120	18.75	20.25		



Laboratory Instrument India