



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
FLARE NUT WRENCHES

**Product Code :**  
LBNY-0029-2200033



**Description :**

Optimum allow-steel structure and head treatment throughout ensure high strength. the long shank of wrench allows extra leverage without compromising safety.  
FINISH: Drop Forged & Chrome Plated.

**Technical Specification :**

SIZE(MM)	SIZE(INCH)	NET WT.IN GMS.	LENGTH(MM)	NET WT. / CTN.(KGS)	GR. WT.(KGS)
8X10	5/16 X 3/8	48 137 24 25.5			
10X11	-- 60	154 24 25.5			
10X12	-- 63	154 25.2 26.7			
11X13	7/16 X 1/2	110 180 44 45.5			
12X14	7/16 X 1/2	110 180 33 34.5			
14X17	-- 154	190 46.2 47.7			
16X17	-- 192	190 28.8 30.3			
17X19	11/16 X 3/4	192 212 28.8 30.3			
18X19	-- 260	212 39 40.5			
19X22	3/4 X 7/8	260 222 26 27.5			
20X22	-- 233	222 23.3 24.8			
21X23	-- 233	222 23.3 24.8			
22X24	-- 350	242 21 22.5			
24X27	-- 428	267 21.4 22.9			
25X28	-- 439	290 17.56 19.06			
30X32	-- 556	291 22.24 23.74			



**Laboratory Instrument India**