



Email : sales@laboratoryinstrumentindia.com

Product Name :
MACHINIST HAMMER WITH FIBRE GLASS HANDLE

Product Code :
LBNY-0030-220009



Description :

MACHINIST HAMMER WITH FIBRE GLASS HANDLE

Technical Specification :

SIZE(GMS)	HANDLE LENGTH(INCHES)	NET. WT. / CTN.(KGS)	GR. WT. / CTN.(KGS)
100	10	20	21.5
200	11	20	21.5
300	12	18	19.5
400	12	20	21.5
500	12	20	21.5
600	14	18	19.5
800	14	20	21.5
1000	14	20	21.5
1500	15	15	16.5
2000	15	20	21.5



Laboratory Instrument India

