



Email : sales@laboratoryinstrumentindia.com

Product Name :
SLEDGE HAMMER WITH FIBER GLASS HANDLE

Product Code :
LBNY-0030-220011



Description :

SLEDGE HAMMER WITH FIBER GLASS HANDLE

Technical Specification :

| SIZE(GMS) | HANDLE LENGTH(INCHES) | NET. WT. / CTN.(KGS) | GR. WT. / CTN.(KGS) |
|-----------|-----------------------|----------------------|---------------------|
| 1 | 10 | 18.2 | 19.7 |
| 1.5 | 10 | 20.475 | 21.975 |
| 2 | 12 | 18.2 | 19.7 |
| 3 | 12 | 20.475 | 21.975 |
| 4 | 12 | 18.2 | 19.7 |
| 6 | 16 | 21.84 | 23.34 |
| 8 | 16 | 21.84 | 23.34 |
| 10 | 16 | 18.2 | 19.7 |
| 12 | 16 | 21.84 | 23.34 |
| 14 | 16 | 19.11 | 20.61 |
| 16 | 16 | 21.84 | 23.34 |
| 18 | 16 | 16.38 | 17.88 |
| 20 | 16 | 18.2 | 19.7 |

Laboratory Instrument India

