



Email : sales@laboratoryinstrumentindia.com

Product Name :
Electrical Tool Box Complete

Product Code :
LBNY-0002-830007



Description :

Electrical Tool Box Complete

Technical Specification :

Electrical tool box complete

This versatile tool kit has screwdrivers and pliers, individually tested for use up to 1000 V ac.

The kit also contains additional extras including sockets, spanners, bit sets, extension bars, hex keys, insulated screwdrivers and VDE pliers in various sizes ensuring that you are covered for a range of applications.

Furthermore, there are useful extras such as scissors, handles and measuring tape to complete the kit, all enclosed with an easy to carry case that is versatile and durable with cut-out slots created for each individual component.



Laboratory Instrument India