



Email : sales@laboratoryinstrumentindia.com

Product Name :
MACHINIST HAMMER HEADS

Product Code :
LBNY-0030-2200019



Description :

Full Drop Forged

Technical Specification :

SIZE(GMS)	NET. WT. / CTN.(KGS)	GR. WT. / CTN.(KGS)
100	20	21.5
200	20	21.5
300	18	19.5
400	20	21.5
500	20	21.5
600	18	19.5
800	20	21.5
1000	20	21.5
1500	15	16.5
2000	20	21.5
3000	18	19.5
4000	16	17.5
5000	20	21.5

Laboratory Instrument India

