



Email : sales@laboratoryinstrumentindia.com

Product Name :
CIRCLIP PLIERS SET 4 PIECE

Product Code :
LBNY-0031-220008



Description :

4pc circlip plier set including internal and external bent and straight galvanized tips, lever handle blue and red dip insulated finger protection.

Features:

Internal and external bent and straight galvanized tips
Insulated finger protection.

Technical Specification :

MATERIAL: Metal
COLOUR: Red
DIMENSIONS: 35 x 20 x 5cm
WEIGHT: 0.83kg.



Laboratory Instrument India